



Recorded: Wednesday, July 31, 2024 at 10:24:33 PM
Heart Rate: 80 BPM
Duration: 30s

Determination: Ectopy
*Kardia Advanced Determination is done on Lead I.

PVC Annotation

Kardia AI can detect the QRS locations within the ECG. It can also identify whether this beat is a PVC or not, the number of PVCs, and the PVC location. This feature is only available for ECGs detected as Sinus Rhythm.

PVC 1

Enhanced Filter, Mains Frequency: 60Hz Scale: 25mm/s, 10mm/mV

