

Patient: XXXXXXXXXX  
Recorded: Tuesday, March 5, 2024 at 11:25:10  
Heart Rate: 129 BPM  
Duration: 30s

**Kardia Advanced** Atrial Fibrillation  
Determination:

Enhanced Filter, Mains Frequency: 60Hz Scale: 25mm/s, 10mm/mV



### R-R Interval Plot

Detection of QRS locations allows Kardia AI to measure the distance between heartbeats or the RR interval. This can be used to review variability of heart rate, which may be useful in understanding heart rate variability, or to visually display irregularity in rate (such as in Afib).

